

 **katena**
DESIGNED FOR SIGHT®

INSTRUMENTS
BIOLOGICS
DEVICES



Knives

Tooke

Corneal Knife
straight blade
with curved cutting edge
K2-3650



Corneal Dissector

For intrastromal dissection.

K2-3660 angled, straight blade
K2-3661 angled, curved blade



K2-3660



K2-3661



Gill

Corneal Knife
curved cutting edge
K2-3700 straight
K2-3710 light curve



Fukasaku

"Hockey" Knife
sharp 3mm front edge
semi-sharp 6mm bottom edge
K2-3735



Spuds • Burrs

Foreign Body Spud

K2-4000 Golf Club, angled, lightly rounded tip

K2-4050 Francis, angled, flag shaped tip

K2-4100 Ellis, lightly curved and rounded tip



K2-4000



K2-4050



K2-4100




Spud & Gouge

in reversible handle

For removing foreign bodies.

K2-4600



Wills Eye Hospital

Foreign Body Spud & Burr
delicate curved spud on one
end and small burr for rust
ring removal on opposite end

K2-4800



Burrs

Algerbrush II

K2-4900 Corneal Rust Ring Remover with protective cap, "AA" battery and two burrs

K2-4903 chuck for K2-4900

K2-4909 burr, 0.5mm


K2-4910 burr, 1.0mm



Pterygium Burr

K2-4912 Pterygium Burr with protective cap, "AA" battery and 3.5mm burr

K2-4913 replacement burr and chuck assembly for K2-4912



Katena

Ophthalmic Burr
for polishing pterygiums

K2-4920 Includes one 3.3mm diameter round diamond burr, two "AA" batteries and storage case.

K2-4923 diamond burr, round 3.3mm diameter

K2-4925 diamond burr, round 5.0mm diameter

K2-4935 diamond burr, disc shaped 5.0mm diameter

K2-4921 optional lithium hi-energy battery, 3.6 volt, for increased torque



Blade Breakers • Diamond Knives

Castroviejo

Blade Breaker and Holder
with sliding lock

K2-5000



Swiss Model

Blade Breaker and Holder
spherical tip

K2-5150




Diamond LRI Step Knife

6-facet blade for bi-directional cutting
single footplate for better visualization
precalibrated depth settings

K2-6518 Precalibrated
depth setting
of 0.6mm

K2-6519 Precalibrated
depth settings
of 0.5mm, 0.55mm
and 0.6mm

K2-6520 Chu LRI Knife,
0.6mm and
2.5mm preset



K2-6519

600 micron knife also available
sterile disposable K20-2560

Diamond Knives

Diamond LRI Micrometer Knife

6-facet blade for bi-directional cutting
 single footplate for better visualization
K2-6522 0 to 1.5mm blade extension
 in 10 micron increments




Diamond AK Knife

double
 vertical and 45°
 cutting edge
K2-6517



Precalibrated depth settings: 0.2, 0.25, 0.3, 0.35, 0.4, 0.6 and 1.5mm



Diamond Micrometer Knife

single
 45° cutting edge
 1mm wide
 200 microns thick
K2-6530



Diamond Micrometer Knife

double
 vertical and 45°
 cutting edge
 200 microns thick
K2-6533



Diamond Micrometer Knife

Thornton "triple"
 vertical and 15°
 cutting edge with
 100 microns base
 cutting edge
K2-6535

Diamond Knives

Blade Gauge

For Katena micrometer knives.

K2-6539



K2-6539

Sterilizing Case

with silicone insert
autoclavable plastic


K2-6555 for one knife

K2-6556 for two knives

K9-2024 for three knives




K2-6555



Diamond Knife

30° cutting edge
1mm wide


K2-6560



Diamond Knife

45° cutting edge
1mm wide


K2-6561



Diamond Knife

tri-facet
1mm wide


K2-6565



Diamond Knife

spear
1mm wide

K2-6567




Diamond Knife


10-facet
1mm wide

K2-6568


Diamond Knives



Diamond Step Knife
single 45° cutting edge
1mm wide
K2-6540




Diamond Step Knife
tri-facet
1mm wide
K2-6546





Diamond Step Knife
10-facet
1mm wide
K2-6548

Precalibrated depth settings: 0.2, 0.3, 0.4, 0.5, 0.6, 0.8 and 4.0mm

Diamond Crescent Knife
crescent shaped blade
2.5mm wide, angled
K2-6575



Diamond Phaco Knife
keratome shaped blade
K2-6578 2.8mm straight
K2-6579 2.8mm angled
K2-6580 3.0mm straight
K2-6581 3.0mm angled
K2-6582 3.2mm straight
K2-6583 3.2mm angled



Diamond Knives • LASEK Trephines

Phaco-4 Diamond Knife

keratome shaped blade
4 cutting edges, angled


K2-6603 2.5mm wide

K2-6605 2.8mm wide

K2-6607 3.0mm wide

K2-6609 3.2mm wide

K2-6611 trapezoid
2.7 to 3.2mm



TRAPEZOID



MICS Diamond Knife


keratome shaped blade
4 cutting edges, angled

K2-6622 trapezoid, 1.8 to 2.2mm wide

K2-6654 trapezoid, 2.1 to 2.4mm wide

K2-6680 1.8mm wide

K2-6682 2.2mm wide



TRAPEZOID



Camellin Style

LASEK Trephine
270° semi-sharp edge
90° cutout for flap hinge

K2-7810 8mm diameter

K2-7812 9mm diameter



Camellin Style

LASEK Trephine &
Alcohol Retaining Well

K2-7830 8mm trephine
and 8.5mm well

K2-7832 9mm trephine
and 9.5mm well



Trabeculectomy Punches

Gass

Corneoscleral Punch
 1.5mm diameter head
 0.75mm deep bite
 rotatable tip for convenient positioning

K2-9500



Luntz-Dodick

Trabeculectomy Punch
 1.0mm diameter head
 with bullet shaped tip
 for scleral tunnel insertion
 0.5mm deep bite
*Cutting head can be preset
 in one of 4 positions to suit
 the surgeon's preference.*

K2-9505 complete

K2-9506 replacement cutting head



Kelly

Descemet's Membrane Punch
 1.0mm diameter head
 0.75mm deep bite

K2-9510



Trabeculectomy Punches

Fukasaku

Micro Trabeculectomy Punch
 0.6mm diameter head
 with bullet shaped tip
 0.3mm x 0.6mm deep bite

K2-9512



Christensen

Punctum Punch
 1mm wide bite
 For opening stenotic punctums.

K2-9525

