

# EASY

refraction unit



- 2 instruments sliding table top
- Power supply built-in
- Simple connection of different types of slit lamps: LED / halogen
- Adjusting knobs for a slit lamp illumination built-in the table top
- Convenient control panel
- Aluminium post with illumination LEDs and universal chart projector adapter
- PA-1000 counterbalanced phoropter arm (equipped with gas spring)
- Steel base powder coated
- Chart projector shelf
- EU Declaration of Conformity provided

# EASY

refraction unit



PA-1000



Control Panel



Auxiliary Outlets Panel



## Dimensions\*

Width: 886 - 1 266mm

Depth: 1 506mm

Height: 1 682mm

Chair height working range: 427 - 627mm

Table top dimensions: 870 x 440mm

Table top height (fixed): 874mm

Table top working range (height adjustable):  
n/a

\*Presented basic parameters may vary depending  
on product version/configuration

## Technical specification\*

Supply voltage: 100-240V AC, 50/60Hz

Power: 280VA

LED lamp power: 9W

Power supply cord: 3m with UniSchuko plug

Max. table top load: 50kg

Max. chair load: 150kg

Max. power of devices supplied from F-type sockets: 450VA

Power supply voltage at the table top sockets (SL/AUX1):  
2.5 - 24V DC

Supply current of the slit lamp fixator:  
0-100mA (LED) or 0-250mA (bulb)

Packaging: 1 pallet, 120x80x160cm, 200kg

# EASY

refraction unit



