

ECO PLUS

refraction unit



- ECO PLUS equipped with electrically elevated COMFORT chair
- 2 instruments sliding table
- Power supply built-in
- Simple connection of different types of slit lamps: LED / halogen
- Convenient control panel and easily accessible accessories drawer
- Aluminium post with illumination LEDs and chart projector shelf
- PA-1000 counterbalanced phoropter arm (equipped with gas spring)
- Steel base powder coated
- EU Declaration of Conformity provided

ECO PLUS

refraction unit




PA-1000



Control Panel



COMFORT chair



Dimensions*

Width: 1 051 - 1 357mm

Depth: 1 414 - 1 522mm

Height: 1 817mm

Chair height working range: 493 - 693mm

Table top dimensions: 880 x 450mm

Table top height (fixed): 878mm

Table top working range (height adjustable):
878 - 978mm

*Presented basic parameters may vary depending
on product version/configuration

Technical specification*

Supply voltage: 100-240V AC, 50/60Hz

Power: 280VA

LED lamp power: 9W

Power supply cord: 3m with UniSchuko plug

Max. table top load: 50kg

Max. chair load: 150kg

Max. power of devices supplied from F-type sockets: 450VA

Power supply voltage at the table top sockets (SL/AUX1):
2.5 - 24V DC

Supply current of the slit lamp fixator:
0-100mA (LED) or 0-250mA (bulb)

Packaging: 1 pallet, 120x80x160cm, ca. 225-270kg (gross)

ECO PLUS

refraction unit

