

SMART PLUS


refraction unit



- SMART PLUS equipped with electrically elevated COMFORT chair
- 2 instruments sliding table top
- Power supply built-in
- Convenient control panel and easily accessible accessories drawer with pivotal arm
- Simple connection of different types of slit lamps: LED / halogen
- Diagonal aluminium post with illumination LEDs and chart projector shelf
- PA-1000 counterbalanced phoropter arm (equipped with gas spring)
- Steel base powder coated
- EU Declaration of Conformity provided

SMART PLUS

refraction unit



Dimensions*

Width: 1 122 - 1 442mm

Depth: 1 536 - 1 644mm

Height: 1 786mm

Chair height working range: 493 - 693mm

Table top dimensions: 900 x 440mm

Table top height (fixed): 878mm

Table top working range (height adjustable):
878 - 978mm

Technical specification*

Supply voltage: 100-240V AC, 50/60Hz

Power: 280VA

LED lamp power: 9W

Power supply cord: 3m with UniSchuko plug

Max. table top load: 50kg

Max. chair load: 150kg

Max. power of devices supplied from F-type sockets: 450VA

Power supply voltage at the table top sockets (SL/AUX1):
2.5 - 24V DC

Supply current of the slit lamp fixator:
0-100mA (LED) or 0-250mA (bulb)

Packaging: 1 pallet, 120x80x160cm, ca. 175-225kg (gross)

*Presented basic parameters may vary depending on product version/configuration

SMART PLUS

refraction unit

