


Laser Indirect Ophthalmoscope (LIO)

Specifications

Product Number	Laser Compatibility	Spot Size*
30903 - H500	OcuLight® GL, GLx, TX, IQ 532™, SL, SLx, IQ 810™	350 µm
13152 - H500	OcuLight SL, SLx, IQ 810	350 µm
13153 - H500	OcuLight SL, SLx, IQ 810	1400 µm
65900 - H500	IQ 577™	350 µm
65515 - H500	IQ 532	350 µm

*on retina using a 20D lens

MicroPulse™ Compatible Lasers: IQ 810, IQ 577, IQ 532, and OcuLight SLx  MicroPulse

1. Shah P, Narendran V, Kalpana N. Large spot transpupillary thermotherapy: A quicker laser for treatment of high risk prethreshold retinopathy of prematurity - a randomized study. Indian J Ophthalmol 2011;59(2):155-8.



Specifications are subject to change without notice. OcuLight, IRIDEX and the IRIDEX logo are registered trademarks and MicroPulse, IQ 810, IQ 577, IQ 532 and TruFocus are trademarks of IRIDEX Corporation. Products are covered by one or more of the following U.S. patents: 5,372,595; 5,511,085; 5,521,932; 5,663,979; 5,982,789; 5,979,554; 6,141,143; 6,144,484; 6,222,869; 6,327,291; 6,377,599; 6,540,391; 6,733,490; 6,800,076; 7,537,593; 7,766,904; 7,771,417; and 7,909,816. Other U.S. and International patents pending. Clinical references available upon request.

 Emergo Europe
Molenstraat 15, 2513 BH The Hague, The Netherlands, Tel.: (31) (0) 70 345-8570, Fax: (31) (0) 70 346-7299

