



# OcuLight<sup>®</sup> SL Laser

## Specifications

<b>Wavelength:</b>	810 nm Infrared
<b>Weight:</b>	14.0 lb (6.4 kg)
<b>Dimensions:</b>	4" H x 12" W x 12" D (10 cm x 30 cm x 30 cm)
<b>Connector Type:</b>	Resistor
<b>Electrical:</b>	100–120 VAC, 50/60 Hz 200–240 VAC, 50/60 Hz
<b>Exposure Duration:</b>	CW-Pulse™: 10–9000 ms
<b>Exposure Interval:</b>	CW-Pulse: 50–1000 ms
<b>Aiming Laser:</b>	Diode laser, 655 nm nominal
<b>Delivery Device</b>	Portable SLA: 0–1300 mW
<b>Power Output:</b>	LIO: 0–1500 mW LIO-LS*: 0–1500 mW EndoProbe <sup>®</sup> : 0–1500 mW G-Probe™: 0–2000 mW DioPexy™: 0–1800 mW OMA: 0–1300 mW

\*Factory installed option



Specifications are subject to change without notice. OcuLight, EndoProbe, IRIDEX and the IRIDEX logo are registered trademarks and G-Probe, CW-Pulse and DioPexy are trademarks of IRIDEX Corporation.

Products are covered by one or more of the following U.S. patents: 5,372,595; 5,511,085; 5,521,932; 5,663,979; 5,982,789; 5,979,554; 6,141,143; 6,144,484; 6,222,869; 6,327,291; 6,377,599; 6,540,391; 6,733,490; 6,800,076; 7,537,593; 7,766,904; 7,771,417; and 7,909,816.

Clinical references available upon request.



Emergo Europe

Molenstraat 15, 2513 BH The Hague, The Netherlands, Tel.: (31) (0) 70 345-8570, Fax: (31) (0) 70 346-7299




# IRIDEX

Elegantly simple solutions™



# A Wide Choice of Delivery Systems

## Specifications

<b>Wavelength:</b>	810 nm Infrared
<b>Weight:</b>	14.0 lb (6.4 kg)
<b>Dimensions:</b>	4" H x 12" W x 12" D (10 cm x 30 cm x 30 cm)
<b>Connector Type:</b>	Resistor
<b>Electrical:</b>	100–120 VAC, 50/60 Hz 200–240 VAC, 50/60 Hz
<b>Exposure Duration:</b>	CW-Pulse™: 10–9000 ms LongPulse™: 10 s–30 min
<b>Exposure Interval:</b>	50–1000 ms
<b>MicroPulse™ Duration:</b>	MicroPulse: 0.1–1.0 ms  MicroPulse
<b>MicroPulse Interval:</b>	MicroPulse: 1.0–10.0 ms
<b>Aiming Laser:</b>	Diode laser, 655 nm nominal
<b>Delivery Device Power Output:</b>	Portable SLA: 0–2000 mW LIO: 0–2000 mW LIO-LS: 0–2000 mW EndoProbe®: 0–2000 mW G-Probe™: 0–3000 mW DioPexy™: 0–2000 mW OMA: 0–2000 mW



Specifications are subject to change without notice. IRIDEX, the IRIDEX logo, OcuLight, and EndoProbe are registered trademarks and CW-Pulse, MicroPulse, LongPulse, G-Probe, and DioPexy are trademarks of IRIDEX Corporation.

Products are covered by one or more of the following U.S. patents: 5,372,595; 5,511,085; 5,521,932; 5,663,979; 5,982,789; 5,979,554; 6,141,143; 6,144,484; 6,222,869; 6,327,291; 6,377,599; 6,540,391; 6,733,490; 6,800,076; 7,537,593; 7,766,904; 7,771,417; and 7,909,816.

 Emergo Europe  
Molenstraat 15, 2513 BH The Hague, The Netherlands, Tel.: (31) (0) 70 345-8570, Fax: (31) (0) 70 346-7299



**Elegantly simple solutions™**

