



Information for Doctors

TRAINING THE BRAIN TO SEE BETTER

RevitalVision™ neural vision therapy
improves communication between the eyes
and the visual cortex.

Proven Results

Improve Vision
**Two Lines on an
Eye Chart**
on Average

**100% Increase
in Contrast Sensitivity**

“I went into this very skeptical, but I did see a positive effect for the patients in the clinical studies. I was impressed that we could improve vision without doing surgery or without changing the optical structure of the eye, but just by improving the brains’ ability to see.”

- Dr. Daniel S. Durrie, MD
Ophthalmologist

For a complete list of clinical trials, and links to medical journals,
please visit, www.RevitalVision.com/Doctors.

RevitalVision Products

REVITALVISION **Post-LASIK** THERAPY

Designed to fine-tune the results of LASIK and other refractive surgeries without corrective lenses or additional surgery. RevitalVision Post-LASIK Therapy has shown to improve visual quality in low light conditions and reduces halos, starbursts and glare while driving at night.

* It is sold only through eye doctors

REVITALVISION **Post-Cataract** THERAPY

Designed to enhance the results of cataract surgery without corrective lenses, additional surgery, drugs or side effects. Patients and clinical trial participants have noticed better vision while driving, less eye fatigue, and have been able to read small print better – such as text on pill bottles or food ingredients.

* It is sold only through eye doctors

REVITALVISION **Vision Improvement** PROGRAM

The Vision Improvement Program can reduce or eliminate the need for reading glasses (presbyopia) and is available online and through eye doctors.

Products Coming Soon:

REVITALVISION **Amblyopia** THERAPY

REVITALVISION **Myopia** THERAPY

REVITALVISION **Pediatric Myopia** THERAPY

REVITALVISION **Sports Vision** PROGRAM

Scientific Background

Based on science created by Nobel Prize winner Dennis Gabor, RevitalVision by NeuroVision is supported by more than two decades of scientific and clinical research.

The technology enhances visual processing in the brain to improve visual acuity and increase contrast sensitivity. Unlike eye exercise programs that focus on relaxing eye muscles, RevitalVision is clinically proven.

During the training sessions the user is presented with a series of precise visual tasks consisting of Gabor patches (Fig 1.0) with subtle differences in orientation, size and contrast. Through repetitive practice the brain is trained to be more efficient and to improve visual processing.

RevitalVision™ algorithms continuously adjust the training sessions based on the patient's performance and neuronal inefficiencies.

More than 3,000 patients in Europe, the U.S. and Asia have been treated.

A more detailed explanation of the science, including links to medical journals and a full list of clinic trails, can be found at www.revitalvision.com/Doctors




Figure 1.0
Gabor patches - Single (top)
and three at 45° angle (bottom)




Figure 2.0
Results published in the
Proceedings of the National
Academy of Sciences,
April 27, 2004

The RevitalVision™ Training System

RevitalVision™ uses neural training to improve vision. The program is:

- Completed at home on a computer
- 2-3 times per week
- Customized to the patient's pace and visual ability
- Professionally monitored by your RevitalVision™ Personal Vision Specialist
- Each of the 20 sessions takes an average of 20 minutes

After they purchase the program from an eye care professional, the product can be shipped directly to the patient's home. Should there be any questions during or after installation, their RevitalVision Personal Vision Specialist will assist them.

RevitalVision monitors progress throughout training, and upon completion the eye doctor will give a final eye exam.

Just like learning to walk, riding a bicycle and learning to swim have long lasting effects, so does RevitalVision. Clinical results show that 85% of what is learned as part of RevitalVision™ therapy is retained.

Safety

Prior to being released for commercial use, reputable eye doctors and eye institutes in the U.S. and across the globe extensively tested RevitalVision by NeuroVision for 20 years. The program is 100% natural and improves vision with no additional surgery, no corrective lenses, no drugs, and no known risks, side effects or complications.

Getting Started

Eye care professionals are encouraged to try the technology before recommending it to customers. Contact RevitalVision to receive your complementary, fully functional, trial version for eye care professionals.

Contact RevitalVision

Phone: (785) 856-0417
Toll Free: 1-866-954-1619
info@revitalvision.com

Founding Ophthalmologists



Dr. Richard L. Lindstrom, MD
Founder, Attending Surgeon of Minnesota Eye Consultants

Dr. Lindstrom is founder and attending surgeon of Minnesota Eye Consultants and Adjunct Professor Emeritus at the University of Minnesota Department of Ophthalmology. He is also medical director of TLC Vision, Sightpath Medical, and Refractec.

He is Chief Medical Editor of the USA and International editions of Ocular Surgery News. He served as President of the American Society of Cataract and Refractive Surgeons.

Dr Lindstrom holds over 30 patents in ophthalmology. He serves on the Board of Directors of AcuFocus, Inc, TLC Vision, Occulogix, eyeonics, Refractec, the Minnesota Medical Foundation, and Inner City Tennis.

Dr. Lindstrom serves on a number of Journal editorial boards and published more than 350 peer reviewed journal articles and 60 book chapters. He is the recipient of numerous awards for distinguished service by national and international ophthalmology associations.



Dr. John D. Hunkeler, MD
Founder, President and CEO of Hunkeler Eye Institute, P.A.

Dr. Hunkeler has more than 30 years experience in cataract and lens implant surgery, corneal and refractive surgery. He is an active investor, advisor and mentor to several start-up and established companies.

Dr. Hunkeler has served as President of the American Society of Cataract and Refractive Surgery, and the Boards of the Outpatient Ophthalmic Surgery Society, and the Eye Bank Association of America.

He has served on several editorial boards professional journals and has numerous scientific publications under his authorship.

Dr. Hunkeler received his Medical Degree from the University of Kansas, School of Medicine and also serves as a Clinical Professor of Ophthalmology at the same location.



Dr. Daniel S. Durrie, MD
Founder, President and CEO of Durrie Vision Center

With more than 30 years experience in refractive and corneal surgery, Kansas City ophthalmologist Dr. Daniel S. Durrie, MD is the leader of the Durrie Vision refractive surgery team. Dr. Durrie has participated in numerous clinical studies and has performed more than 40,000 refractive surgery procedures.

Dr. Durrie's scientific articles and papers have appeared in numerous professional medical publications. He has served on numerous editorial boards for several professional journals.

Dr. Daniel Durrie received his medical doctorate and completed his residency at the University of Nebraska. He serves as Clinical Professor and Director of Refractive Surgery Services at the University of Kansas Medical Center.

RevitalVision is:

Clinically Proven and
Tested Worldwide

Published in
Peer Reviewed Medical Journals

Supported by a
World Renowned Medical Advisory Board



RevitalVision, LLC
1617 St. Andrews, Suite 210
Lawrence, Kansas 66047

Phone: (785) 856-0417
Toll Free: 1-866-954-1619
info@RevitalVision.com

www.revitalvision.com