


GoVu™


PACHYMETER



Hand-held and easy to use

Instantly and accurately measures  
central corneal thickness

Battery-operated with  
charging cradle included



Integrated 20MHz transducer for excellent accuracy

Full scan control available at fingertips

Ergonomic shape comfortably fits either hand

Easy navigation of scan results and features

Large screen for viewing results and battery status

## Charge. Then Go Vu Central Corneal Thickness.

### Specifications

Length	7.5"
Weight	3.6 oz
Display	1.1" Color LCD
Battery	Rechargeable Lithium-Ion 3.7V / 750 mAh
Charger	Input: 90-264 VAC (Universal) Output: 5 VDC / 1A
Transducer	20 MHz frequency
Measurement Range	300 - 1000 $\mu$ m
Clinical Accuracy	Within 5 $\mu$ m
Measurement Clock	160 MHz
Averaging	25 measurements per displayed result
Precision	1 $\mu$ m



Sonomed Escalon™

Sonomed Escalon  
1979 Marcus Avenue C105  
Lake Success, NY 11042 USA

800-227-1285 or 516-354-0900  
info@sonomedinc.com  
www.sonomedescalon.com